# 

### PRODUCT

## PRINCIPLE Sterile Nitrile Gloves

#### DESCRIPTION

Integrity Cleanroom PRINCIPLE Sterile Nitrile Gloves are hand specific, with a longer length beaded cuff which has been designed for extra protection. The textured surface and enhanced grip makes these sterile nitrile gloves ideal for handling delicate items.

Our gloves are powder and latex free, assisting in reducing the risk of an allergic reaction.

Sterilized by minimum 2.5Mrad Gamma radiation.

#### FEATURES

- Sterile
- Class 100/ISO 5
- 100% nitrile butadiene rubber
- Powder and latex free
- Length: 305mm / 12 inch
- Longer length for enhanced protection
- Textured surface for easy handling

PRODUCT CODE	DESCRIPTION	SIZE	QUANTITY (per case)
602-0400	PRINCIPLE Sterile Nitrile Glove	6	1 pair/pouch   10 pairs/outer   200 pairs/case
602-0401	PRINCIPLE Sterile Nitrile Glove	6.5	1 pair/pouch   10 pairs/outer   200 pairs/case
602-0402	PRINCIPLE Sterile Nitrile Glove	7	1 pair/pouch   10 pairs/outer   200 pairs/case
602-0403	PRINCIPLE Sterile Nitrile Glove	7.5	1 pair/pouch   10 pairs/outer   200 pairs/case
602-0404	PRINCIPLE Sterile Nitrile Glove	8	1 pair/pouch   10 pairs/outer   200 pairs/case
602-0405	PRINCIPLE Sterile Nitrile Glove	8.5	1 pair/pouch   10 pairs/outer   200 pairs/case
602-0406	PRINCIPLE Sterile Nitrile Glove	9	1 pair/pouch   10 pairs/outer   200 pairs/case
602-0407	PRINCIPLE Sterile Nitrile Glove	10	1 pair/pouch   10 pairs/outer   200 pairs/case

To request a quotation or for more information, please call **+1 512-713-0127** email **sales@integritycleanroom.com** or visit **www.integritycleanroom.com** 

IMPORTANT: This data sheet and its contents (the "Information") belong to Integrity Cleanroom or are licensed to it. No licence is granted for the use of it other than for information purposes in connection with the products to which it relates. No licence of any intellectual property rights is granted. The Information is subject to change without notice and replaces all data sheets previously supplied. The Information supplied is believed to be accurate but Integrity Cleanroom assumes no responsibility for its accuracy or completeness, any error in or omission from it or for any use made of it. Users of this data sheet should check for themselves the Information and the subtability of their purpose and not make any assumptions based on information included or omitted. Liability for loss or damage resulting from negligence or where Integrity Cleanroom was aware of the possibility of such loss or damage arising) is excluded. © 2025 Integrity Cleanroom.



#### SPECIFICATIONS

PROPERTIES	RESULT			
Surface	Palm textured			
Color	White			
Length	305mm / 12 inch			
Thickness	5.0 mil (±1)			
Tensile strength	Min 15 Mpa			
Elongation	600%			
Packing style	1 pair per pouch, 10 pairs per pack, 20 packs per case (200 pairs total)			
Physical Properties	Meets the requirements of the Glove Standard EN420:2003+A1:2009 and the European Medical Gloves Standard EN455 (Parts 1, 2, 3, 4) & EN374-:2016 (Part 1, 2, 3, 4 and 5). Product fulfills the cleanroom specifications of nitrile gloves as per ISO 23464:2020			
AQL	AQL 1.5			
Material	100% Nitrile Butadiene Rubber			
Shelf life	Five (5) years from date of irradiation			
Storage instructions	Store in a cool dry place (< 104°F) away from direct sunlight and fluorescent light			
Batch Traceability	Individual bags of pairs marked with batch numbers			
Liquid Particle Counts (as determined by IES-RP-C005.4)	< 1200 particle counts cumulative at 0.5micron			
Organics (IES-RP-CC005.4)	No Silicone Oil / DOP / Amides detected by FTIR Spectroscopy			
Non-Volatile Residue	<10 microgram/ sq. cm			
ESD Parameters (ASTM D-257)	Surface Resistance - ≤10E10 Ohms/sq. Decay Time (Sec) - ≤0.5 Tribo Charge (Volts) - ≤20			
Dimension Criteria	Length: 305 ±10 mm Palm Width: 77 ±5 mm to 120 ±5 mm Single Wall thickness Finger thickness: 0.18 ± 0.02 mm Palm thickness: 0.12 ± 0.02 mm			
	Palm thickness: 0.12 ± 0.02 mm Cuff thickness: 0.09 ± 0.02 mm			

Ionic Burden (microgram/ sq cm)(as determined by IES-RP-C005.4)						
Fluoride Chloride Phosphate Sodium	<0.030 <1.200 <0.800 <0.050	Sulphate Nitrate Bromide Magnesium Potassium	<0.800 <1.000 <0.050 <0.004 <0.050			

#### To request a quotation or for more information, please call **+1 512-713-0127** email **sales@integritycleanroom.com** or visit **www.integritycleanroom.com**

IMPORTANT: This data sheet and its contents (the "Information") belong to Integrity Cleanroom or are licensed to it. No licence is granted for the use of it other than for information purposes in connection with the products to which it relates. No licence of any intellectual property rights is granted. The Information is subject to change without notice and replaces all data sheets previously supplied. The Information supplied is believed to be accurate but Integrity Cleanroom assumes no responsibility for its accuracy or completeness, any error in or omission from it or for any use made of it. Users of this data sheet should check for themselves the Information and the suitability of their purpose and not make any assumptions based on information included or omitted. Liability for loss or damage resulting from any reliance on the Information or use of it (including liability resulting from negligence or where Integrity Cleanroom was aware of the possibility of such loss or damage arising) is excluded. © 2025 Integrity Cleanroom.